

OASIS<sup>®</sup>



---

PRODUCT CATALOG

*International*

A decorative border at the bottom of the page consisting of a series of parallel diagonal lines in a dark gray color.





## ABOUT OASIS

OASIS Medical, Inc., is dedicated to equipping eye care providers with cost-effective solutions to improve patient management and practice growth. OASIS supplies high-quality surgical and dry eye products to eye care professionals that provide solutions to chronic conditions and make a meaningful difference in patients' quality of life.

### **A company founded on exceptional standards**

Founded in 1987 by Norman Delgado and headquartered near Los Angeles, OASIS Medical, Inc. is a specialty manufacturer of eye care products and surgical supplies, instruments, and disposables. For 25 years, OASIS has been manufacturing a broad range of products for use in cataract, refractive, dry eye, and retinal disease for domestic and international sale.

### **An expanding portfolio**

The founding vision of OASIS was to ensure that the company continually addresses the changing needs of surgeons, clinicians, and patients. OASIS offers a portfolio of exceptional products, while also improving the way in which products are packaged, delivered, or sold. The original product portfolio of OASIS comprised single-use surgical items. Since then, we have expanded to offer a wide range of more than 300 ophthalmic and optometric accessories, supplies and instruments.

### **Looking forward**

OASIS efficiently delivers solution-focused products to eye care practices. OASIS will continue to provide innovative, differentiated products for eye care providers and consumers.



## OUR CORPORATE MESSAGE

### **Our Philosophy**

OASIS is a continuous improvement company deeply committed to truth in its business approach, relationships, and products based on a foundation of accountability. We welcome feedback from customers and our team to evolve as a proactive company.

### **Our Investment**


OASIS invests in the long term through our strategic partnerships with eye care providers. This investment helps the doctor improve the quality of the patient's eye care experience. As consultants we are willing to invest time and equip the doctor with how our approach can assist their practice. This strategy allows OASIS to thrive by being relevant through focusing on the patient care experience and efficiency of the individual practice.

### **Our Solution Based Products**

Our values of integrity and quality guide our commitment to solution-based products that enhance patient care. From the beginning, OASIS has been a practice-based company focused on building dry eye and surgical practices through the various tools we offer.

### **Our Approach**

OASIS is committed to offering informed, relevant solutions based on truth. To do this, we listen first and ask questions to understand the issues. Our professionalism is demonstrated by our ability to think differently. We lead with a consultative approach that is flexible and seeks to add value. It looks to address issues rather than only sell a product. It is strategic rather than tactical. Our goal is to create a partnership between OASIS and doctors so that both parties win.



# DRY EYE PRODUCTS

The OASIS Dry Eye product line contains punctum and intracanalicular plugs, instruments, diagnostics, and consumer products.

## EXTRA INFORMATION

### Patient Education

OASIS offers patient-friendly educational materials about common ocular conditions, disease, and our products. In support of eye care practices, we offer patient education aids, instructions for using OASIS products, patient consent forms, publications regarding OASIS products, and information about reimbursement. Our goal is to assist eye care providers in improving patient outcomes and augmenting practice growth.

For more information:  
Visit the RESOURCES tab at  
[www.oasismedical.com](http://www.oasismedical.com)

## TABLE OF CONTENTS

<a href="#">  FORM FIT Intracanalicular Plugs</a>	06
<a href="#">  SOFT PLUG Silicone Punctum Plugs</a>	07
<a href="#">  SOFT PLUG Flow Control Plugs</a>	08
<a href="#">  SOFT PLUG Collagen Plugs</a>	08
<a href="#">  SOFT PLUG Extended Duration Plugs</a>	09
<a href="#">  Dry Eye Marketing Materials</a>	10
<a href="#">  Instruments &amp; Diagnostics</a>	10



# SURGICAL PRODUCTS

OASIS offers an expansive line of reusable and disposable surgical products across cataract, refractive, retinal, and lacrimal surgery applications.



## TABLE OF CONTENTS

OASIS Iris Expander	11
Premier Edge Scalpels	12 - 15
Disposable Instruments	16 - 26
Collagen Shield	27
Viscoelastic	28
LASIK	29 - 31
Accessory Products	32 - 33
Product Index	34 - 35
International Distribution	36
Shipping Policy	37
Ordering Information	38
Quality Assurance Information	39

## EXTRA INFORMATION

OASIS is diligent in offering high quality, cost-effective products that can improve surgical efficiencies, ensuring high standards of safety, performance and compliance.



# FORM FIT®

## INTRACANALICULAR PLUGS

### ONE SIZE. ONE SOLUTION.




Form Fit is truly a one size plug solution. Once inserted into the canaliculus, Form Fit hydrates over a 10 minute period. As it hydrates, Form Fit increases in size until it completely fills the vertical canalicular cavity.

### HOW IT WORKS

Form Fit is made of a hydrogel material that expands into a soft, pliable, gelatinous material when it contacts tear film. It takes approximately 10 minutes for the plug to become completely hydrated. Form Fit securely conforms to the punctal cavity so there is no need for a cap or anchor system on the end of the plug. This greatly increases patient comfort and decreases foreign body sensation. Form Fit does not begin to swell until it makes contact with tear film reducing the potential for product waste. Form Fit is made from a hydrogel material that is safe and biocompatible.

### PRODUCT FEATURES


- Hydrates over a period of 10 minutes
- Requires no special storage
- Easily removed via irrigation of the tear duct
- Biocompatible, no foreign body sensation
- Conveniently pre-loaded on inserter
- 0.3mm, one size fits all



Product Name	Description	Ref #
FORM FIT Intracanalicular Plug 10pk	One size	6303
FORM FIT Intracanalicular Plug 2pk	One size	6303-T







## SOFT PLUG® SILICONE PUNCTUM PLUGS

The distinctive low-profile cap combined with a softer, more flexible silicone maximizes patient comfort. The inserter provides greater control and eliminates accidental release of the plug. Soft Plug silicone plugs are removable.

Available preloaded in 2-pack and 6-pack dispenser box or non-loaded Economy Pack.

Product Name	Description	Size	Ref #
SOFT PLUG Preloaded Silicone 2pk	Micro	0.4mm	6610
SOFT PLUG Preloaded Silicone 2pk	Mini	0.5mm	6611
SOFT PLUG Preloaded Silicone 2pk	Petite	0.6mm	6612
SOFT PLUG Preloaded Silicone 2pk	Small	0.7mm	6613
SOFT PLUG Preloaded Silicone 2pk	Medium	0.8mm	6614

Product Name	Description	Size	Ref #
SOFT PLUG Preloaded Silicone 6pk	Micro	0.4mm	6610-D
SOFT PLUG Preloaded Silicone 6pk	Mini	0.5mm	6611-D
SOFT PLUG Preloaded Silicone 6pk	Petite	0.6mm	6612-D
SOFT PLUG Preloaded Silicone 6pk	Small	0.7mm	6613-D
SOFT PLUG Preloaded Silicone 6pk	Medium	0.8mm	6614-D

Product Name	Description	Size	Ref #
SOFT PLUG Silicone Economy Pack	Micro	0.4mm	6610-EP
SOFT PLUG Silicone Economy Pack	Mini	0.5mm	6611-EP
SOFT PLUG Silicone Economy Pack	Petite	0.6mm	6612-EP
SOFT PLUG Silicone Economy Pack	Small	0.7mm	6613-EP
SOFT PLUG Silicone Economy Pack	Medium	0.8mm	6614-EP

- Economy pack includes 2 inserters and 10 trays of 2 plugs each. Trays equipped with an alignment slot for quick loading.



## SOFT PLUG® FLOW CONTROL PLUGS

2 plugs per box ■

Preloaded on inserter ■

Silicone ■

For partial occlusion with limited tear drainage ■




Product Name	Description	Size	Ref #
SOFT PLUG Flow Control 2pk	Petite	0.6mm	6612-FC
SOFT PLUG Flow Control 2pk	Small	0.7mm	6613-FC
SOFT PLUG Flow Control 2pk	Medium	0.8mm	6614-FC

## SOFT PLUG® COLLAGEN INTRACANALICULAR PLUGS

The collagen plugs have a smooth proprietary finish that simplifies insertion of the plug. Collagen plugs are ideal for short-term, diagnostic, and postsurgical occlusion. Plugs are 2mm in length and absorbable within 2-5 days.

60 plugs per box.



Product Name	Size	Ref #
SOFT PLUG Collagen Plugs	0.2mm	6502
SOFT PLUG Collagen Plugs	0.3mm	6503
SOFT PLUG Collagen Plugs	0.4mm	6504

**THIS PRODUCT IS NOT SOLD IN ALL COUNTRIES.  
PLEASE CONTACT OASIS FOR MORE INFORMATION.**







## SOFT PLUG® EXTENDED DURATION PLUGS

Extended Duration absorbable plugs are designed to provide patient comfort up to 3 months. This longlasting plug is effective and economical for transient dry eye symptoms after surgery. Extended Duration plugs are 2mm in length and sterile.

Product Name	Size	Ref #
SOFT PLUG Extended Duration Plugs 20pk	0.2mm	6402-A
SOFT PLUG Extended Duration Plugs 20pk	0.3mm	6403-A
SOFT PLUG Extended Duration Plugs 20pk	0.4mm	6404-A
SOFT PLUG Extended Duration Plugs 20pk	0.5mm	6405-A

Product Name	Size	Ref #
SOFT PLUG Extended Duration Plugs 40pk	0.2mm	6402
SOFT PLUG Extended Duration Plugs 40pk	0.3mm	6403
SOFT PLUG Extended Duration Plugs 40pk	0.4mm	6404
SOFT PLUG Extended Duration Plugs 40pk	0.5mm	6405

Product Name	Size	Ref #
SOFT PLUG Extended Duration Plugs Variety Pack	1 pair 0.2mm, 3 pairs 0.3mm, 4 pairs 0.4mm, and 2 pairs 0.5mm	6400V



## DRY EYE MARKETING MATERIALS

Product Name	Ref #
Dry Eye Patient Brochure (English)	6520
Dry Eye Patient Brochure (Spanish)	6520S
Patient Education Video	6530
Punctal Occlusion Training DVD	6532
Dry Eye & Surgery Patient Brochure	6540
Oasis Wall Poster	6542-I
Dry Eye Symptoms Wall Poster	6543

OASIS offers educational materials in the form of brochures, displays, posters, and videos for use in your practice and to distribute to your patients.

To view and download digital versions of all resources listed, click on the RESOURCES tab at [www.oasismedical.com](http://www.oasismedical.com)

## INSTRUMENTS & DIAGNOSTICS


Product Name	Description	Ref #
Sizing Gauge	.5mm and .6mm tips	6625
Sizing Gauge	.7mm and .8mm tips	6627
Sizing Gauge Set	Includes both instruments & sterilization tray	6630




Product Name	Description	Ref #
Straight Jeweler's Forceps	Straight #3	6643
Curved Jeweler's Forceps	Curved #7	6647
Grooved Forceps	For plug insertion	6648
Castroviejo Lacrimal Dilator	Double Ended, 14cm	6640



Product Name	Description	Ref #
Lacrimal Irrigation Kit	One 3ml syringe and one 25g x 1 1/8" lacrimal cannula, curved, open end. 10 per box.	4003-IK





## OASIS IRIS EXPANDER

### EXPAND YOUR VIEW. SUSTAIN YOUR ACCESS.

The OASIS Iris Expander is a 7mm diameter polypropylene ring that expands the iris and maintains access and visibility throughout the surgical procedure.

The disposable device is indicated in cases where Miosis and Intraoperative Floppy Iris Syndrome are present.

Each OASIS Iris Expander is easy to insert and remove and gentle to fragile iris tissue. Each OASIS Iris Expander is packaged sterile with its own easily loaded disposable injector.

### PRODUCT FEATURES

- 7mm diameter
- Contains ring and disposable injector
- Polypropylene material
- Sterile
- Single use
- Easy to insert and remove

Product Name	Ref #
OASIS Iris Expander 5pk	9700
OASIS Iris Expander 1pk	9700S





## PREMIER EDGE® SCALPELS

6 per box ■

Disposable plastic handle ■

Stainless steel micro blades ■


Sterile blister pack ■



Angle	Ref #
15°	PE3015-SL

### PREMIER EDGE SAFETY INCISION SCALPEL

Sterile, 6 per box.



Width	Ref #
2.0mm	PE3720A-SL

### PREMIER EDGE SAFETY CRESCENT SCALPEL

Sterile, 6 per box.



Width	Ref #
2.2mm	PE4822-SL
2.4mm	PE4824-SL
2.65mm	PE4826-SL
2.75mm	PE4827-SL
2.85mm	PE4828-SL
3.0mm	PE4830-SL
3.2mm	PE4832-SL

### PREMIER EDGE SAFETY THIN LINE SLIT SCALPEL

Sterile, 6 per box.





### PREMIER EDGE SAFETY MVR SCALPEL

Sterile, 6 per box.

Gauge	Ref #
20g	PE3520-SL



### PREMIER EDGE SIDE PORT SCALPEL

Sterile, 6 per box.

Width	Ref #
1.0mm	PE3010
1.2mm	PE3012
1.4mm	PE3014



### PREMIER EDGE RESTRICTED DEPTH STAB SCALPEL

Sterile, 6 per box.

Depth	Ref #
3.0mm	PE3015-3
5.0mm	PE3015-5



### PREMIER EDGE INCISION/STAB SCALPEL

Sterile, 6 per box.

Angle	Ref #
15°	PE3015
22.5°	PE3022
30°	PE3030
45°	PE3045



### PREMIER EDGE IMPLANT SCALPEL

Sterile, 6 per box.

Width	Ref #
4.0mm	PE3440
5.2mm	PE3452



Angle	Width	Ref #
Angled	1.4mm	PE3814
Angled	1.7mm	PE3817
Angled	2.2mm	PE3822
Angled	2.4mm	PE3824
Angled	2.5mm	PE3825
Angled	2.65mm	PE3826
Angled	2.75mm	PE3827
Angled	2.85mm	PE3828
Angled	3.0mm	PE3830
Angled	3.2mm	PE3832
Straight	3.2mm	PE3832S

Width	Ref #
2.2mm	PE4822-TL
2.4mm	PE4824-TL
2.65mm	PE4826-TL
2.75mm	PE4827-TL
2.85mm	PE4828-TL
3.0mm	PE4830-TL
3.2mm	PE4832-TL

Width	Ref #
2.65mm	PE3226
2.75mm	PE3227
2.85mm	PE3228
3.0mm	PE3230

## PREMIER EDGE SLIT SCALPEL

Sterile, 6 per box.



## PREMIER EDGE THIN LINE SLIT SCALPEL

Sterile, 6 per box.



## PREMIER EDGE CLEAR CORNEA PARALLEL SIDE SCALPEL

Sterile, 6 per box.







### PREMIER EDGE CLEAR CORNEA TRAPEZOID SCALPEL

Sterile, 6 per box.

Width	Ref #
1.8 - 2.2mm	PE3301
2.5 - 3.0mm	PE3302



### PREMIER EDGE CLEAR CORNEA SCALPEL

Sterile, 6 per box.

Width	Ref #
2.2mm	PE3122
2.5mm	PE3125
2.65mm	PE3126
2.75mm	PE3127
2.85mm	PE3128
3.2mm	PE3132



### PREMIER EDGE MVR SCALPEL

Sterile, 6 per box.

Gauge	Angle	Ref #
19g	Straight	PE3519
20g	Straight	PE3520
23g	Straight	PE3523
25g	Straight	PE3525
19g	Angled	PE3519A
20g	Angled	PE3520A
23g	Angled	PE3523A



### PREMIER EDGE TUNNEL SCALPEL

Sterile, 6 per box.

Type	Width/Angle	Ref #
Round	2.2mm/Angled	PE3622
Crescent	2.0mm/Angled	PE3720A
Crescent	2.0mm/Straight	PE3720



## DISPOSABLE INSTRUMENTS

OASIS supplies an extensive line of cannulas, needles, and cystotomes. We also offer customized delivery options in order to provide cost savings to practices.




## ANESTHESIA

### RETROBULBAR ANESTHESIA NEEDLE

Sterile, 10 per box.


Description	Ref #
25g (0.50) x 1-1/2" (38mm)	4004



### RETROBULBAR ATKINSON ANESTHESIA NEEDLE

Sterile, 10 per box.


Description	Ref #
23g (0.60) x 1-1/2" (38mm)	4004C
25g (0.50) x 1-1/2" (38mm)	4004B



### PERIBULBAR ATKINSON ANESTHESIA NEEDLE


Sterile, 10 per box.

Description	Ref #
25g (0.50) x 7/8" (22mm)	4004D
23g (0.60) x 1-1/4" (32mm)	4004F
25g (0.50) x 1-1/4" (32mm)	4004M



The Atkinson shortened single bevel decreases possibility of retrobulbar hemorrhage.






## SUB TENON CANNULA

Sterile, 10 per box.

Description	Ref #
19g (1.10) x 1" (25mm) 25mm radius, flattened to .66mm	4005




## TRI-PORT SUB TENON CANNULA

Sterile, 10 per box.

Description	Ref #
19g (1.10)x 1" (25mm) Two 0.4mm opposite side ports	4005TP


## IRRIGATION/AIR INJECTION &amp; VISCOELASTIC




## ANTERIOR CHAMBER CANNULA

Sterile, 10 per box.

Description	Ref #
25g (0.50) x 7/8" (22mm) 4mm bend	4001AF
27g (0.40) x 7/8" (22mm) 4mm bend	4002F
30g (0.30) x 7/8" (22mm) 4mm bend	4001F



Description	Ref #
30g (0.30) x 7/8" (22mm) 6mm bend	4001
25g (0.50) x 7/8" (22mm) 6mm bend	4001A
27g (0.40) x 7/8" (22mm) 6mm bend	4002



## CURVED CANNULA

Sterile, 10 per box.

Description	Ref #
25g (0.50) x 1-1/8" (29mm)	4003






## ANTERIOR CHAMBER WASHOUT CANNULA

Sterile, 10 per box.


Description	Ref #
19g (1.10) x 1" (25mm) 8mm bend	4036




## ANTERIOR CHAMBER/VISCOELASTIC CANNULA

Sterile, 10 per box.

Description	Ref #
23g (0.60) x 7/8" (22mm) 10mm bend	4001B
27g (0.40) x 7/8" (22mm) 11mm bend	4002A




Description	Ref #
23g (0.60) x 7/8" (22mm) 8mm bend, Thin Wall	4001CTW
27g (0.40) x 7/8" (22mm) 8mm bend	4002B



## VISCOELASTIC CANNULA

Sterile, 10 per box.

Description	Ref #
27g (0.40) x 7/8" (22mm) 9mm bend, 45° Angle	4002L





## HYDRODISSECTION

### HORIZONTAL-FLOW HYDRODISSECTION CANNULA

Sterile, 10 per box.

Description	Ref #
25g (0.50) x 7/8" (22mm) 8mm bend, Flattened tip	4036A
27g (0.40) x 7/8" (22mm) 8mm bend, Flattened tip	4036B






VERTICAL-FLOW HYDRODISSECTION CANNULA

Sterile, 10 per box.


Description	Ref #
25g (0.50) x 7/8" (22mm) 4mm bend, Flattened tip	4036C
27g (0.40) x 7/8" (22mm) 4mm bend, Flattened tip	4036D



HELSINKI HYDRODISSECTION CANNULA

Sterile, 10 per box. Curved with flattened tip. (Horizontal Flow)

Description	Ref #
25g (0.50) x 7/8" (22mm)	4036F
27g (0.40) x 7/8" (22mm)	4036G




CHANG HYDRODISSECTION CANNULA

Sterile, 10 per box.

90° angled flattened tip provides a broad stream of fluid.

The point of the beveled tip is used to engage and rotate the nucleus to ensure complete separation.

Description	Ref #
27g (0.40) x 3/4" (19mm) 1mm bend	4036J



SAUTER HYDRODISSECTION CANNULA

Sterile, 10 per box. Curved with flattened tip (Vertical Flow)


Description	Ref #
25g (0.50) x 7/8" (22mm)	4036M
27g (0.40) x 7/8" (22mm)	4036N



## CORTICAL CLEAVING HYDRODISSECTOR

Sterile, 10 per box.


Description	Ref #
27g (0.40) x 7/8" (22mm) 11mm bend, flat end	4036Q
27g (0.40) x 7/8" (22mm) 11mm bend, round end	4036P



## U-SHAPED HYDRODISSECTION CANNULA

Sterile, 10 per box.


Description	Ref #
25g (0.50) x 7/8" (22mm) 2.5mm length x 2.0mm wide hook	4016AH



## KELLAN HYDRODISSECTION CANNULA

Sterile, 10 per box. Curved, 20mm radius. Flattened 3mm from end.

Description	Ref #
25g (0.50) x 7/8" (22mm)	4068





## CAPSULOTOMY


## FORMED IRRIGATING CYSTOTOME

Sterile, 10 per box.

Description	Ref #
25g (0.50) x 7/8" (22mm)	4017
27g (0.40) x 7/8" (22mm)	4017B
25g (0.50) x 5/8" (16mm)	4017G
27g (0.40) x 5/8" (16mm)	4017IO
30g (0.30) x 5/8" (16mm)	4017I
25g (0.50) x 5/8" (16mm) Radius of curvature 1/2"	4017IE








## STRAIGHT IRRIGATING CYSTOTOME

Sterile, 10 per box.


Description	Ref #
25g (0.50) x 7/8" (22mm)	4017A
27g (0.40) x 5/8" (16mm)	4017IS



## REVERSE FORMED IRRIGATING CYSTOTOME

Sterile, 10 per box.


Description	Ref #
25g (0.50) x 5/8" (16mm)	4017L



## ANGLED IRRIGATING CYSTOTOME

Sterile, 10 per box.

Description	Ref #
25g (0.50) x 5/8" (16mm) 12mm bend	4017M
27g (0.40) x 5/8" (16mm) 12mm bend	4017N



## CAPSULORHEXIS CYSTOTOME

Sterile, 10 per box. Formed.

Designed with cutting edge from front to back (push-pull) for circular tear capsulotomies. Tip is angled 90° and then rotated 90°.


Description	Ref #
25g (0.40) x 5/8" (16mm)	4017S
27g (0.40) x 5/8" (16mm)	4017R



## F-SL CURVE IRRIGATING CYSTOTOME

Sterile, 10 per box.

Description	Ref #
25g (0.40) x 5/8" (16mm) 19mm radius	4017Z




## CORTEX EXTRACTORS

## U-SHAPED CORTEX EXTRACTOR

Sterile, 10 per box.


Description	Ref #
21g (0.80) x 7/8" (22mm) 2.0mm L x 2.25mm W hook	4016
23g (0.60) x 7/8" (22mm) 2.0mm L x 2.25mm W hook	4016F
27g (0.40) x 7/8" (22mm) 2.5mm L x 2.0mm W hook	4016B



## SIMCOE CORTEX EXTRACTOR

Sterile, 10 per box. Curved, 0.3mm top port


Description	Ref #
23g (0.60) x 1/2" (13mm)	4008



## CHARLEUX CORTEX EXTRACTOR

Sterile, 10 per box. Curved, 10mm radius, 45° tip


Description	Ref #
23g (0.60) x 1/2" (13mm)	4026




## MICRO J-SHAPED CORTEX EXTRACTOR

Sterile, 10 per box. 1.5mm length x 1.75mm width hook.

Description	Ref #
25g (0.50) x 7/8" (22mm)	4016G





## CAPSULE POLISHERS



## KRATZ CAPSULE POLISHER

Sterile, 10 per box. Roughened tip.


Description	Ref #
21g (0.80) x 7/8" (22mm) 9mm bend, 0.3mm top port	4027
23g (0.60) x 7/8" (22mm) 9mm bend, 0.3mm top port	4027A
21g (0.60) x 7/8" (22mm) 6mm bend, open end	4027B



## TERRY CAPSULE POLISHER

Sterile, 10 per box. 10mm silicone sleeve.

Description	Ref #
23g (0.60) x 7/8" (22mm) 7mm bend	4028B
27g (0.40) x 7/8" (22mm) 6mm bend	4028



## OLIVE-TIP ANGLED CAPSULE POLISHER

Sterile, 10 per box. 2mm olive tip.

Description	Ref #
25g (0.50) x 1" (25mm) 10mm bend	4027E




## IRRIGATING & ASPIRATING CANNULAS

### SIMCOE DOUBLE I & A CANNULA

Sterile

Description	Ref #
23g/23g (0.60mm/0.60mm) 0.30mm aspiration port	4007




## LENS LOOPS

### IRRIGATING VECTIS


For use with gravity fluid source

Description	Ref #
24g (0.55) x 5.5mm (Three 0.25mm ports)	4029A



## INFUSION CANNULA

Description	Ref #
2.0mm	4042-2S
4.0mm	4042-4S



## ANTERIOR CHAMBER MAINTAINER

20g (0.90) ■

8.5" silicone tubing with 3mm tip ■


Provides irrigation for chamber maintenance during anterior surgery ■



Product Name	Gauge	Ref #
Anterior Chamber Maintainer	20g (0.90)	4300P




## LACRIMAL



## CURVED OPEN END LACRIMAL CANNULA

Sterile, 10 per box.


Description	Ref #
21g (0.80) x 1-1/2" (38mm)	6005



## CURVED CLOSED END LACRIMAL CANNULA

Sterile, 10 per box.


Description	Ref #
25g (0.50) x 1-1/2" (38mm) 0.3mm port facing down 1mm from tip	6007



## STRAIGHT CLOSED END LACRIMAL CANNULA

Sterile, 10 per box.

Description	Ref #
21g (0.80) x 1-1/2" (38mm) 0.3mm port 1mm from tip	6006



## STRAIGHT STEPPED LACRIMAL CANNULA

Sterile, 10 per box.


Description	Ref #
19g (1.10) x 1-1/2" (38mm) 12.5mm tip stepped down to 23g (0.60)	6008






## LACRIMAL INTUBATION SET W/ RETRIEVAL DEVICE

Description	Ref #
Probes 23g (0.60) x 4-1/4"	
Medical grade silicone tubing length 14"	
.025" O.D. / .012" I.D.	6004
Stainless steel malleable wire	
Probe end has special notch for retrieval instrument to catch	





## LACRIMAL INTUBATION SET

Description	Ref #
Probes 23g (0.60) x 7"	
Medical grade silicone tubing length 11-3/4"	6001
.025" O.D. / .012" I.D.	



## DCR SET

Description	Ref #
Straight	
Medical grade silicone tubing length 11-3/4" Diameter 0.03"	6002
Curved	
Medical grade silicone tubing length 11-3/4" Diameter 0.03"	6003
45° angle 15mm from distal end of probes	



THIS PRODUCT IS NOT SOLD IN ALL COUNTRIES.  
PLEASE CONTACT OASIS FOR MORE INFORMATION.





## SOFT SHIELD® COLLAGEN CORNEAL SHIELD

- Widest range of dissolution rates available
- Thin, purified collagen
- 14.5mm diameter shield
- Accurate and consistent dissolution times

Provides ocular comfort, protection and lubrication following surgery, traumatic or non-traumatic corneal conditions.

Product Name	Duration	Ref #
SOFT SHIELD Quick Shield Collagen Corneal Shield	Less than 12 hours	7002
SOFT SHIELD Collagen Corneal Shield	12 hours non-patched, 24 hours patched	7012
SOFT SHIELD Collagen Corneal Shield	24 hours non-patched, 72 hours patched	7024
SOFT SHIELD Collagen Corneal Shield	72 hours non-patched	7072

**THIS PRODUCT IS NOT SOLD IN ALL COUNTRIES.  
PLEASE CONTACT OASIS FOR MORE INFORMATION.**



## VISCO SHIELD® VISCOELASTIC

Dispenser box, 6 syringes per box ■

2ml syringe with 23g x 7/8" thin wall cannula ■

HPMC (hydroxypropyl methylcellulose) ■

Stores at room temperature ■



VISCO SHIELD is a dispersive HPMC viscoelastic composed of a sterile, high molecular weight, hydroxypropyl methylcellulose solution for use in anterior segment surgical procedures including cataract extraction and intraocular lens implantation.

Product Name	Description	Ref #
VISCO SHIELD® Viscoelastic	20,000 cps	5122
VISCO SHIELD® Viscoelastic	8,000 cps	51082

**THIS PRODUCT IS NOT SOLD IN ALL COUNTRIES.  
PLEASE CONTACT OASIS FOR MORE INFORMATION.**

## VISCO SHIELD® TOPICAL

6 per box ■

1.5ml filled syringe with 23g x 7/8" thin wall cannula ■

HPMC (hydroxypropyl methylcellulose) ■

Preservative-free, stores at room temperature ■




Designed for topical use to hydrate and protect the outer surface of the eye, Visco Shield Topical is preservative-free and designed to keep the cornea hydrated for up to 20 minutes. Visco Shield Topical adheres to the cornea and forms a thin, optically neutral clear film.

Product Name	Description	Ref #
VISCO SHIELD Topical	8,000 cps	5200




# LASIK CANNULAS



## LIN LASIK CANNULA

Sterile, 10 per box.


Description	Ref #
25g x 7/8", 9mm flattened angle, two 0.3mm ports	4036L



## BURATTO LASIK CANNULA

Sterile, 10 per box.


Description	Ref #
25g, 11mm angle, one front and two side ports	4036LB



## LINDSTROM CANNULA

Sterile, 10 per box.


Description	Ref #
25g, 45° angle, 5/8" radius, three 0.3mm ports	4036LR



## SLADE LASIK IRRIGATOR

Sterile, 10 per box.

Description	Ref #
26g x 1", 9mm formed extension with spatulated tip	4036LS



## LASIK OPEN END CANNULA


Sterile, 10 per box.

Description	Ref #
25g x 7/8", Formed 7mm from Open End Six 0.15mm side ports at 2mm, 4mm & 6mm from end	4225



## SOFT CELL® PVA FOAM

Fiber free	■
Highly absorbent	■
Fast wicking	■
Enhanced fluid capacity	■



## Description

## Ref #

High density spears designed specifically for LASIK. Fiber-free, Fast Wicking.

0525

SOFT CELL  
SURGICAL SPEARS

Sterile, 150 per box



## Description

## Ref #

Designed for use in delicate surgical procedures

0529

SOFT CELL  
WIPE + WICK

1 sterile Wipe + Wick per package,  
20 packages per box.



## Description

## Ref #

Eye Pledget


0527

SOFT CELL  
EYE PLEDGET



Description	Ref #
9mm	0528

## SOFT CELL LASIK SHIELD



Description	Ref #
(Chayet) 10mm internal diameter	0522-5
(Chayet/Katzen) 11.5mm internal diameter. Allows visualization of axis marks at 3:00/9:00 for torsional registration. Reduces tear film debris in the interface. Reduces reflux pooling at the hinge during ablation.	0523-5

## SOFT CELL LASIK EYE DRAIN

Sterile, 20 per box



## COLVARD PUPILLOMETER

**Accurately measure scotopic pupil size and corneal diameter of cataract and refractive candidates prior to surgery.**

Utilizes light amplification technology to allow visualization of the patient's pupils in a darkened exam room. Includes storage case, lithium battery and instruction manual.



Product Name	Description	Ref #
Colvard Pupilometer	Light weight, hand held, battery operated instrument.	0401A

## DISC FLUID FILTER

Filters out particulate and biological contaminate matter larger than 0.22  $\mu$ .

Bidirectional filter: female luer lock and male luer slip with screw-in lock to secure cannula in place.

For use between syringe and needle.



Product Name	Description	Ref #
Disc Fluid Filter	50 per box	0601

## OASIS EYE CAUTERY

Fine tip ■

Battery operated ■

Low temperature ■

Range: 1100°F  $\pm$  50°F ■



Product Name	Description	Ref #
OASIS Eye Cautery	10 per box	3501-A





## IRIS RETRACTORS

- 6 pouches containing 5 hooks each per box
- To temporarily stretch and retain the iris during surgery
- Flexible blue polypropylene hooks for better visibility
- Flat silicone stopper adjusts easily to secure the hook

### Product Name

### Ref #

Iris Retractors

4465



## STERILIZATION TRAY

- Tray for day sterilization
- Holds 5 scalpels
- 2 bars

### Product Name

### Ref #

Sterilization Tray

0701



## PRE-INKED CORNEAL MARKING PEN

- Gentian Violet pre-inked pen
- Fine tip
- 10 pens per box

### Product Name

### Ref #

Corneal Marking Pen

0442



# PRODUCT INDEX

Ref #	Product Name	Page
0401A	Colvard Pupillometer	32
0442	Pre-Inked Corneal Marking Pen	33
0525	SOFT CELL Surgical Spear	30
0527	Eye Pledget	30
0528	LASIK Shield	31
0529	SOFT CELL Wipe + Wick	30
0522-5 - 0523-5	SOFT CELL LASIK Eye Drain	31
0601	Disc Fluid Filter	32
0701	Sterilization Tray	33
PE3015-SL	Premier Edge Safety Incision	12
PE3720A-SL	Premier Edge Safety Crescent	12
PE4822-SL - PE4832-SL	Premier Edge Safety Thin Line Slit	12
PE3520-SL	Premier Edge Safety MVR	13
PE3010 - PE3014	Premier Edge Side Port	13
PE3015 - PE3045	Premier Edge Stab/Incision	13
PE3015-3 - PE3015-5	Premier Edge Restricted Depth Stab	13
PE3122 - PE3132	Premier Edge Clear Cornea	15
PE3226 - PE3230	Premier Edge Clear Cornea Parallel Side	14
PE3301 - PE3302	Premier Edge Clear Cornea Trapezoid	15
PE3440 - PE3452	Premier Edge Implant	13
PE3519 - PE3525	Premier Edge MVR	15
PE3622 - PE3720A	Premier Edge Tunnel	15
PE3814 - PE3832S	Premier Edge Slit	14
PE4822-TL - PE4832-TL	Premier Edge Thin Line Slit	14
3501-A	OASIS Eye Cautery	32
4001AF - 4002	Anterior Chamber Cannula	17

Ref #	Product Name	Page
4001CTW	Thin Wall Anterior Chamber Cannula	18
4001B - 4002B	Anterior Chamber/Viscoelastic Cannula	18
4002L	Viscoelastic Cannula	18
4003	Curved Cannula	17
4003-IK	Lacrimal Irrigation Kit	10
4004	Retrobulbar Anesthesia Needle	16
4004B - 4004C	Retrobulbar Atkinson Anesthesia Needle	16
4004D - 4004M	Peribulbar Atkinson Anesthesia Needle	16
4005TP	Tri-Port Sub Tenon Cannula	17
4005	Sub Tenon Cannula	17
4007	Simcoe Double I & A Cannula	24
4008	Simcoe Cortex Extractor	22
4016 - 4016F	U-Shaped Cortex Extractor	22
4016AH	U-Shaped Hydrodissection Cannula	20
4016G	Micro J-Shaped Cortex Extractor	22
4017A - 4017IS	Straight Irrigating Cystotome	21
4017 - 4017IE	Formed Irrigating Cystotome	20
4017L	Reverse Formed Irrigating Cystotome	21
4017M - 4017N	Angled Irrigating Cystotome	21
4017R - 4017S	Capsulorhexis Cystotome	21
4017Z	F-SL Curve Irrigating Cystotome	22
4026	Charleux Cortex Extractor	22
4027 - 4027B	Kratz Capsule Polisher	23
4027E	Olive-Tip Angled Capsule Polisher	23
4028 - 4028B	Terry Capsule Polisher	23
4029A	Lens Loop Irrigating Vectis	24



# PRODUCT INDEX

Ref #	Product Name	Page
4036	Anterior Chamber Washout Cannula	18
4036A - 4036B	Horizontal-Flow Hydrodissection Cannula	18
4036C - 4036D	Vertical-Flow Hydrodissection Cannula	19
4036F - 4036G	Helsinki Hydrodissection Cannula	19
4036J	Chang Hydrodissection Cannula	19
4036L	Lin LASIK Cannula	29
4036LB	Buratto LASIK Cannula	29
4036LR	Lindstrom LASIK Irrigator	29
4036LS	Slade LASIK Irrigator	29
4036M - 4036N	Sauter Hydrodissection Cannula	19
4036P - 4036Q	Cortical Cleaving Hydrodissector	20
4042-2S - 4042-4S	Infusion Cannula	24
4068	Kellan Hydrodissection Cannula	20
4225	LASIK Open End Cannula	29
4300P	Anterior Chamber Maintainer	24
4465	Iris Retractors	33
5122 - 51082	VISCO SHIELD Viscoelastic	28
5200	VISCO SHIELD Topical	28
6001	Lacrimal Intubation Set	26
6002 - 6003	DCR Sets	26
6004	Lacrimal Intubation Set w/ Retrieval Device	26
6005	Curved Open End Lacrimal Cannula	25
6006	Curved Closed End Lacrimal Cannula	25
6007	Straight Closed End Lacrimal Cannula	25
6008	Straight Stepped Lacrimal Cannula	25
6303-T	FORM FIT Plugs 2-pk	6

Ref #	Product Name	Page
6303	FORM FIT Plugs 10-pk	6
6400V	SOFT PLUG Ext Dur Plugs Variety Pack	9
6402-A - 6405-A	SOFT PLUG Ext Dur Plugs 20pk	9
6402 - 6405	SOFT PLUG Ext Dur Plugs 40pk	9
6502 - 6504	SOFT PLUG Collagen Plugs	8
6520 - 6544S	Dry Eye Marketing Materials	10
6610 - 6614	SOFT PLUG Silicone 2pk	7
6610-D - 6614-D	SOFT PLUG Silicone 6pk	7
6610-EP - 6614-EP	SOFT PLUG Silicone Econ Pack	7
6612-FC - 6614-FC	SOFT PLUG Flow Control	8
6625 - 6630	Sizing Gauges	10
6640	Castroviejo Lacrimal Dilator	10
6643	Straight Jeweler's Forceps	10
6647	Curved Jeweler's Forceps	10
6648	Grooved Forceps	10
7002 - 7072	SOFT SHIELD Corneal Shield	27
9700 - 9700S	OASIS Iris Expander	11





## INTERNATIONAL DISTRIBUTION

Customers outside the USA, please order from the OASIS distributor in your country.

For more information regarding your distributor, please e-mail: [international@oasismedical.com](mailto:international@oasismedical.com)

**If you are interested in becoming an OASIS distributor, please contact the International Department: [international@oasismedical.com](mailto:international@oasismedical.com)**

OASIS currently has distributors in most countries worldwide:

Algeria	Croatia	Hungary	Netherlands	Slovakia
Argentina	Cyprus	Iceland	New Zealand	Slovenia
Australia	Czech Republic	India	Norway	South Africa
Austria	Denmark	Indonesia	Pakistan	Spain
Belgium	Dominican Republic	Israel	Peru	Sweden
Bolivia	Ecuador	Italy	Philippines	Switzerland
Brazil	Egypt	Japan	Poland	Taiwan
Bulgaria	El Salvador	Jordan	Portugal	Thailand
Canada	Finland	Korea	Puerto Rico	Turkey
Caribbean	France	Kuwait	Romania	UAE
Chile	Germany	Lebanon	Russia	Ukraine
China	Greece	Lithuania	Saudi Arabia	United Kingdom
Colombia	Guatemala	Malaysia	Serbia	Venezuela
Costa Rica	Hong Kong	Mexico	Singapore	Vietnam



**Donawa Lifescience Consulting Srl**

Piazza Albania 10

00153 Rome, Italy



## SHIPPING POLICY

International orders are normally shipped via FedEx to include shipping charges. If a preferred carrier other than FedEx is requested, the distributor will be asked to provide the contact details of the carrier and account number. Customers in countries outside the U.S. are charged separately for all VAT, duties, and brokerage fees. Next-day service to Puerto Rico is available but not guaranteed. We cannot be responsible for delays in transit.

In most situations, in-stock product orders received before 2pm PST (UTC/GMT -8 hours) Monday-Friday, excluding holidays, will ship on the next business day. Exceptions are when backorders occur or when customers are on a prepayment basis. Prepayment orders will ship once payment is received. International orders will receive a pro forma confirmation of product ordered on the following business day. The pro forma will include shipping charges. Once the order has shipped, tracking numbers will be sent via email.

If a distributor requests bundled shipment from various purchase orders, a five (5) day window will be placed on products withheld from inventory on any of those orders. New distributors may require an eight (8) week window on shipment of first orders.

## TRADEMARKS AND COPYRIGHTS

**Trademarks owned and/or registered to OASIS Medical, Inc. are:**

OASIS®, SOFT PLUG®, SOFT SHIELD®, FORM FIT®, SOFT CELL®, VISCO SHIELD®, PREMIER EDGE®, PREMIER SHIELD®, Oasis TEARS®, Oasis TEARS® PLUS, Oasis TEARS® MULTIDOSE

This catalog is copyrighted and may not be reproduced without the consent of OASIS Medical Inc.



## ORDERING INFORMATION

Orders can be placed 24 hours via our website, fax, or the phone message center.

### CUSTOMER SERVICE HOURS

**7:00 AM - 4:30 PM** (PST),  
Monday - Friday, Excluding holidays.

**PAYMENT TERMS:** Net 30 days

**SHIPPING:** FOB Glendora, California  
Insurance and shipping charges added to invoice.

**USA:** 800.528.9786  
800.631.7210 FAX

**WORLDWIDE:** 909.305.5400  
909.305.6897 FAX

**ONLINE:** [www.oasismedical.com](http://www.oasismedical.com)

## RETURN GOODS POLICY

Unopened items returned for credit must be accompanied by an RGA (return goods authorization number) and the reason for the return. An RGA can be requested through customer service.

### Allowance Schedule for acceptable returns

Sellable merchandise returned in factory condition within 60 days	Full Credit
Sellable merchandise returned in factory condition within 60-90 days	Cost less 40%
Sellable merchandise returned in factory condition 90 days to 1 year	Cost less 60%
Merchandise returned one year after invoice date	No Credit
Merchandise unsuitable for resale	No Credit
Merchandise defective at time of receipt by customer	Full Credit
Merchandise shipped in error	Full Credit

### Merchandise that will not be accepted:

- "Special Order" or "Custom Products"
- Merchandise not in standard OASIS packaging
- Opened products
- Discontinued products

**NOTE:** Customer must prepay all freight charges on stock returns except for defective material or merchandise shipped in error.



# QUALITY ASSURANCE WORLDWIDE

OASIS is registered with the FDA as a:

- Medical Device Establishment
- Drug Establishment

The OASIS Quality Management System is ISO 13485-2003 certified by DQS Medizinprodukte GmbH.



## CE Marked Product Lines

Silicone Punctum Plugs	Lacrimal Instruments
Cannulas, Needles, Cystotomes	Astigmatome
HPMC Viscoelastic	Microkeratome Blades
PVA Foam Sponge Products	Extended Duration Plug
Pupillometer	Corneal Marking Instrument
Syringe Disc Filters	Bimanual Irrigation/Aspiration Sets
Eye Cautery	Microsurgical Knives
Hydrogel Plug	

## STERILIZATION METHOD

Products are sterilized by ethylene oxide, gamma radiation, electron beam, or moist heat according to EN 556 (European) guidelines and ANSI/AAMI/ISO 11135, 11137 or 17665 standards for medical device sterilization as applicable, to a sterility assurance level (SAL) of  $10^{-6}$ .

## PRODUCT LABELING

OASIS uses internationally recognized symbols on its product labeling.

**Product name and description** — **SOFT CELL®**  
PVA FOAM SURGICAL SPEAR  
LANCETTE CHIRURGICALE OPHTHALMIQUE EN MOUSSE PVA  
CHIRURGISCHER SPEER AUS PVA-SCHAUMSTOFF  
HEMOSTETA QUIRÚRGICA DE ESPUMA PVA  
MICROSPUGNETTE TRIANGOLARI AD ALTO ASSORBIMENTO SU BASTONCINO

**European Union CE Mark and notified body ID number** — **CE 0297**

**Method of Sterilization:**  
**R** = gamma or e-beam radiation  
**EO** = ethylene oxide  
 = steam or dry heat

**Do not reuse** — **SINGLE USE**

**Label part number and country of origin** — **P627 7/98 MADE IN USA**

**Reference Catalog #** — **REF [Part No.] 0525**

**Lot or batch code** — **LOT F0698C**

**Use by** — **2001-07**

**Expiration Date - year/month** — **514 SO. VERMONT AVE. GLENDORA, CA 91741**

**SN: Serial Number** — **Quantity: 10 spears**

**Attention, see instructions for use** — **Do not use if package is opened or damaged.**

**OASIS®**







514 S. Vermont Ave. Glendora, California 91741

[international@oasismedical.com](mailto:international@oasismedical.com) | (909) 305-5400 | [www.oasismedical.com](http://www.oasismedical.com)